

## Product datasheet for **GA214750**

### Mouse Trarg1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                               |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                  |
| Description:  | Trarg1 CRISPRa kit - CRISPR gene activation of mouse trafficking regulator of GLUT4 (SLC2A4)<br>1                                                                                                                                                                                                             |
| Symbol:       | Trarg1                                                                                                                                                                                                                                                                                                        |
| Components:   | GA214750G1, Trarg1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214750G2, Trarg1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA214750G3, Trarg1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q8C838</a>                                                                                                                                                                                                                                                                                        |
| Synonyms:     | C130069F04Rik; DSPB1; Lost1                                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                       |

#### Product images:



[View online »](#)