

Product datasheet for **GA214773**

Mouse Wdcp activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdcp CRISPRa kit - CRISPR gene activation of mouse WD repeat and coiled coil containing
Symbol:	Wdcp
Components:	GA214773G1, Wdcp gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214773G2, Wdcp gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214773G3, Wdcp gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6NV72</u>
Synonyms:	A430092C04
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)