

Product datasheet for **GA214776**

Mouse *Kcns3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcns3</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, delayed-rectifier, subfamily S, member 3
Symbol:	<i>Kcns3</i>
UniProt ID:	Q8BQZ8
Synonyms:	D12Erttd137e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)