

Product datasheet for **GA214785**

Mouse *Kcnh5* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnh5</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, subfamily H (eag-related), member 5
Symbol:	<i>Kcnh5</i>
UniProt ID:	Q920E3
Synonyms:	<i>Eag2</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)