

Product datasheet for **GA214794**

Mouse *Serpina3f* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Serpina3f CRISPRa kit – CRISPR gene activation of mouse serine (or cysteine) peptidase inhibitor, clade A, member 3F
Symbol:	Serpina3f
Components:	GA214794G1, Serpina3f gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214794G2, Serpina3f gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214794G3, Serpina3f gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q80X76
Synonyms:	2A1; BC049975
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

