

Product datasheet for **GA214811**

Mouse Zfp367 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp367 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 367
Symbol:	Zfp367
Components:	GA214811G1, Zfp367 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214811G2, Zfp367 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214811G3, Zfp367 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q0VDT2</u>
Synonyms:	8030486J21Rik; ZFF29a; ZFF29b; Znf367
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)