


Product datasheet for **GA214825**

Mouse **Lrit1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Lrit1 CRISPRa kit - CRISPR gene activation of mouse leucine-rich repeat, immunoglobulin-like and transmembrane domains 1
Symbol:	Lrit1
Components:	GA214825G1, Lrit1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214825G2, Lrit1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214825G3, Lrit1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8K099
Synonyms:	BC032270; Lrrc21
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)