

## Product datasheet for **GA214839**

### Mouse Kctd12 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Kctd12 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 12 |
| Symbol:       | Kctd12  |
| UniProt ID:   | <a href="#">Q6WVG3</a>  |
| Synonyms:     | AU046135; AW538430; Pfet1; Pfetin   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                     |

#### Product images:



[View online >](#)