

Product datasheet for **GA214875**

Mouse Lamp3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lamp3 CRISPRa kit - CRISPR gene activation of mouse lysosomal-associated membrane protein 3
Symbol:	Lamp3
UniProt ID:	Q7TST5
Synonyms:	1200002D17Rik; Cd208; DC-LAMP; DCLAMP; LAMP; TSC403
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)