

Product datasheet for **GA214935**

Mouse Lipk activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Lipk CRISPRa kit - CRISPR gene activation of mouse lipase, family member K |
| Symbol: | Lipk |
| UniProt ID: | Q8BM14 |
| Synonyms: | 9930115F20; BC055815; Lip12 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)