


## Product datasheet for **GA214970**

### Mouse Plekhm3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Plekhm3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family M, member 3
Symbol:	Plekhm3
Components:	GA214970G1, Plekhm3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA214970G2, Plekhm3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA214970G3, Plekhm3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q8BM47</a></u>
Synonyms:	9430067K14Rik; A230102O09Rik; AI449570; DAPR; Plekhm11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)