

Product datasheet for GA214974

Mouse *Asic4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Asic4</i> CRISPRa kit - CRISPR gene activation of mouse acid-sensing (proton-gated) ion channel family member 4
Symbol:	<i>Asic4</i>
UniProt ID:	Q7TNS7
Synonyms:	<i>Accn4</i> ; <i>BNAC4</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)