


Product datasheet for **GA215004**

Mouse **Ldlrad3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Ldlrad3 CRISPRa kit - CRISPR gene activation of mouse low density lipoprotein receptor class A domain containing 3
Symbol:	Ldlrad3
Components:	GA215004G1, Ldlrad3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215004G2, Ldlrad3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215004G3, Ldlrad3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2AR95
Synonyms:	6430500P08; AI194318; Lrad3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)