

Product datasheet for **GA215006**

Mouse Wdr76 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr76 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 76
Symbol:	Wdr76
Components:	GA215006G1, Wdr76 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215006G2, Wdr76 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215006G3, Wdr76 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6PWY4
Synonyms:	5830411K18Rik; 6030457D09; AA517836
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)