

Product datasheet for **GA215025**

Mouse Pabpc4l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pabpc4l CRISPRa kit - CRISPR gene activation of mouse poly(A) binding protein, cytoplasmic 4-like
Symbol:	Pabpc4l
Components:	GA215025G1, Pabpc4l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215025G2, Pabpc4l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215025G3, Pabpc4l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	C330050A14Rik; EG241989
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)