

Product datasheet for **GA215035**

Mouse Dnai3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr63 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 63
Symbol:	Dnai3
Components:	GA215035G1, Dnai3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215035G2, Dnai3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215035G3, Dnai3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RY71</u>
Synonyms:	4931433A13Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

