

Product datasheet for **GA215036**

Mouse *Lrrc7* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Lrrc7</i> CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 7
Symbol:	<i>Lrrc7</i>
Components:	GA215036G1, <i>Lrrc7</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215036G2, <i>Lrrc7</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215036G3, <i>Lrrc7</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80TE7</u>
Synonyms:	AK122498; densin; mKIAA1365
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)