


Product datasheet for **GA215057**

Mouse Kank4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kank4 CRISPRa kit - CRISPR gene activation of mouse KN motif and ankyrin repeat domains 4
Symbol:	Kank4
Components:	GA215057G1, Kank4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215057G2, Kank4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215057G3, Kank4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6P9J5
Synonyms:	Ankrd38; AU021071; BC060737; C130031J23
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)