


Product datasheet for **GA215060**

Mouse Slc35d1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc35d1 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 35 (UDP-glucuronic acid/UDP-N-acetylgalactosamine dual transporter), member D1
Symbol:	Slc35d1
Components:	GA215060G1, Slc35d1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215060G2, Slc35d1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215060G3, Slc35d1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2AKQ0</u>
Synonyms:	AI834976; C330011J09; mKIAA0260; UGTREL7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)