

Product datasheet for **GA215079**

Mouse *Lrrc38* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lrrc38 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 38
Symbol:	Lrrc38
Components:	GA215079G1, Lrrc38 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215079G2, Lrrc38 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215079G3, Lrrc38 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A2VDH3
Synonyms:	A230053A07Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)