

Product datasheet for **GA215096**

Mouse Kctd8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kctd8 CRISPRa kit - CRISPR gene activation of mouse potassium channel tetramerisation domain containing 8
Symbol:	Kctd8
Components:	GA215096G1, Kctd8 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215096G2, Kctd8 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215096G3, Kctd8 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q50H33
Synonyms:	A730087N02Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)