

## Product datasheet for **GA215097**

### Mouse Tecrl activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tecrl CRISPRa kit - CRISPR gene activation of mouse trans-2,3-enoyl-CoA reductase-like
Symbol:	Tecrl
Components:	GA215097G1, Tecrl gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215097G2, Tecrl gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215097G3, Tecrl gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<a href="#">Q8BFZ1</a>
Synonyms:	D330017N19Rik; Srd5a2l2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)