


Product datasheet for **GA215098**

Mouse Tmprss11f activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmprss11f CRISPRa kit - CRISPR gene activation of mouse transmembrane protease, serine 11f
Symbol:	Tmprss11f
Components:	GA215098G1, Tmprss11f gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215098G2, Tmprss11f gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215098G3, Tmprss11f gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BHM9</u>
Synonyms:	4631408C20; 4732406D01Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)