

Product datasheet for GA215099

Mouse Ugt2b35 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ugt2b35 CRISPRa kit – CRISPR gene activation of mouse UDP glucuronosyltransferase 2 family, polypeptide B35
Symbol:	Ugt2b35
Components:	GA215099G1, Ugt2b35 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215099G2, Ugt2b35 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215099G3, Ugt2b35 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	C730031G17
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

