

Product datasheet for **GA215134**

Mouse Cfap100 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cfap100 CRISPRa kit - CRISPR gene activation of mouse cilia and flagella associated protein 100
Symbol:	Cfap100
Components:	GA215134G1, Cfap100 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215134G2, Cfap100 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215134G3, Cfap100 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80VN0</u>
Synonyms:	AW048498; C030041G11Rik; C230069K22Rik; Ccdc37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)