

Product datasheet for **GA215144**

Mouse *Plxna4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Plxna4 CRISPRa kit - CRISPR gene activation of mouse plexin A4
Symbol:	Plxna4
UniProt ID:	<u>Q80UG2</u>
Synonyms:	9330117B14; mKIAA1550; Plxa4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)