

Product datasheet for **GA215167**

Mouse Siglecg activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Siglecg CRISPRa kit - CRISPR gene activation of mouse sialic acid binding Ig-like lectin G
Symbol:	Siglecg
UniProt ID:	<u>Q80ZE3</u>
Synonyms:	9830164H23; A630096C01Rik; mSiglec-G; Siglec-G; Siglec10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)