

Product datasheet for **GA215168**

Mouse Shank1 activation kit by CRISPRa

Product data:

Product Type: CRISPRa Kits

Description: Shank1 CRISPRa kit - CRISPR gene activation of mouse SH3 and multiple ankyrin repeat domains 1


Symbol: Shank1

Components: GA215168G1, Shank1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)
GA215168G2, Shank1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)
GA215168G3, Shank1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)
1 CRISPRa-Enhancer vector, SKU GE100056
1 CRISPRa scramble vector, SKU GE100077

UniProt ID: D3YZU1

Format: 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)