

Product datasheet for **GA215169**

Mouse Zfp473 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp473 CRISPRa kit - CRISPR gene activation of mouse zinc finger protein 473
Symbol:	Zfp473
Components:	GA215169G1, Zfp473 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215169G2, Zfp473 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215169G3, Zfp473 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BI67</u>
Synonyms:	D030014N22Rik; zfp-100; Znf473
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)