

Product datasheet for GA215200

Mouse Wdr17 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr17 CRISPRa kit – CRISPR gene activation of mouse WD repeat domain 17
Symbol:	Wdr17
Components:	GA215200G1, Wdr17 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215200G2, Wdr17 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215200G3, Wdr17 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	3010002I12Rik; B230207L18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

