

Product datasheet for **GA215210**

Mouse Pkd1l3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pkd1l3 CRISPRa kit – CRISPR gene activation of mouse polycystic kidney disease 1 like 3
Symbol:	Pkd1l3
Components:	GA215210G1, Pkd1l3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215210G2, Pkd1l3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215210G3, Pkd1l3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q2EG98</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

