

Product datasheet for **GA215232**

Mouse Tnfaip8l3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnfaip8l3 CRISPRa kit - CRISPR gene activation of mouse tumor necrosis factor, alpha-induced protein 8-like 3
Symbol:	Tnfaip8l3
Components:	GA215232G1, Tnfaip8l3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215232G2, Tnfaip8l3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215232G3, Tnfaip8l3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3TBL6</u>
Synonyms:	9930029P06Rik; Gm517
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)