

Product datasheet for GA215291

Mouse Pwwp3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pwwp3b CRISPRa kit – CRISPR gene activation of mouse PWWP domain containing 3B
Symbol:	Pwwp3b
Components:	GA215291G1, Pwwp3b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215291G2, Pwwp3b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215291G3, Pwwp3b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q4VA55</u>
Synonyms:	9430086K2IRik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

