

Product datasheet for **GA215300**

Mouse Tceanc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tceanc CRISPRa kit - CRISPR gene activation of mouse transcription elongation factor A (SII) N-terminal and central domain containing
Symbol:	Tceanc
Components:	GA215300G1, Tceanc gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215300G2, Tceanc gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215300G3, Tceanc gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3US16
Synonyms:	A630018P17Rik; Gm390
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)