

Product datasheet for GA215304

Mouse Ssh3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ssh3 CRISPRa kit - CRISPR gene activation of mouse slingshot protein phosphatase 3
Symbol:	Ssh3
UniProt ID:	Q8K330
Synonyms:	BC028922; mSSH-3L; SSH-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)