

Product datasheet for **GA215306**

Mouse Spag4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Spag4 CRISPRa kit - CRISPR gene activation of mouse sperm associated antigen 4
Symbol:	Spag4
Components:	GA215306G1, Spag4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215306G2, Spag4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215306G3, Spag4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9JIF2
Synonyms:	1700041K21Rik; mKIAA4118; MNCb-0953; Sun4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)