

Product datasheet for **GA215308**

Mouse Pcmt2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcmt2 CRISPRa kit - CRISPR gene activation of mouse protein-L-isoaspartate (D-aspartate) O-methyltransferase domain containing 2
Symbol:	Pcmt2
Components:	GA215308G1, Pcmt2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215308G2, Pcmt2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215308G3, Pcmt2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BHD8</u>
Synonyms:	5330414D10Rik; AI314201
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)