


Product datasheet for **GA215316**

Mouse Slc36a2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc36a2 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 36 (proton/amino acid symporter), member 2
Symbol:	Slc36a2
Components:	GA215316G1, Slc36a2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215316G2, Slc36a2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215316G3, Slc36a2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8BHK3
Synonyms:	A530067G19Rik; PAT2; Tramd1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)