

Product datasheet for GA215327

Mouse Kcne2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Kcne2 CRISPRa kit - CRISPR gene activation of mouse potassium voltage-gated channel, Isk-related subfamily, gene 2
Symbol:	Kcne2
UniProt ID:	Q9D808
Synonyms:	2200002I16Rik; AW048273; MiRP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)