

Product datasheet for **GA215355**

Mouse Cd300lf activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cd300lf CRISPRa kit - CRISPR gene activation of mouse CD300 molecule like family member F
Symbol:	Cd300lf
UniProt ID:	Q6S1Q7
Synonyms:	CLIM1; CLM-1; CLM1; Digr2; F730004D16Rik; IgSF13; IREM1; LMIR3; Pigr3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)