

Product datasheet for **GA215360**

Mouse Trpv3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trpv3 CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily V, member 3
Symbol:	Trpv3
UniProt ID:	<u>Q8K424</u>
Synonyms:	1110036I10Rik; A1644701; Nh; VRL3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)