

Product datasheet for **GA215367**

Mouse **Dusp15** activation kit by **CRISPRa**

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Dusp15 CRISPRa kit - CRISPR gene activation of mouse dual specificity phosphatase-like 15 |
| Symbol: | Dusp15 |
| UniProt ID: | Q8R4V2 |
| Synonyms: | AI851682; LMW-DSP10; T-DSP10 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)