

Product datasheet for **GA215373**

Mouse Ap1s3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ap1s3 CRISPRa kit - CRISPR gene activation of mouse adaptor-related protein complex AP-1, sigma 3
Symbol:	Ap1s3
UniProt ID:	Q7TN05
Synonyms:	1190009B22; Jr2; [s]3A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)