

## Product datasheet for **GA215404**

### Mouse Olfr735 activation kit by CRISPRa

#### Product data:

|                      |                                                                                                                                                                                                                                                                                                                  |
|----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type:</b> | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| <b>Description:</b>  | Olfr735 CRISPRa kit – CRISPR gene activation of mouse olfactory receptor 735                                                                                                                                                                                                                                     |
| <b>Symbol:</b>       | Olfr735                                                                                                                                                                                                                                                                                                          |
| <b>Components:</b>   | GA215404G1, Olfr735 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215404G2, Olfr735 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215404G3, Olfr735 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| <b>Synonyms:</b>     | MOR243-1                                                                                                                                                                                                                                                                                                         |
| <b>Format:</b>       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:

