


## Product datasheet for **GA215480**

### Mouse Olfr122 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Olfr122 CRISPRa kit - CRISPR gene activation of mouse olfactory receptor 122                                                                                                                                                                                                                                     |
| Symbol:       | Olfr122                                                                                                                                                                                                                                                                                                          |
| Components:   | GA215480G1, Olfr122 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215480G2, Olfr122 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215480G3, Olfr122 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | MOR263-10                                                                                                                                                                                                                                                                                                        |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)