


## Product datasheet for **GA215580**

### Mouse Olfr318 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Olfr318 CRISPRa kit - CRISPR gene activation of mouse olfactory receptor 318                                                                                                                                                                                                                                     |
| Symbol:       | Olfr318                                                                                                                                                                                                                                                                                                          |
| Components:   | GA215580G1, Olfr318 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215580G2, Olfr318 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215580G3, Olfr318 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | MOR285-1; RP23-128D9.10                                                                                                                                                                                                                                                                                          |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)