

## Product datasheet for **GA215618**

### Mouse Olfr1020 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                     |
|---------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                        |
| Description:  | Olfr1020 CRISPRa kit - CRISPR gene activation of mouse olfactory receptor 1020                                                                                                                                                                                                                                      |
| Symbol:       | Olfr1020                                                                                                                                                                                                                                                                                                            |
| Components:   | GA215618G1, Olfr1020 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215618G2, Olfr1020 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215618G3, Olfr1020 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8VFK7</a></u>                                                                                                                                                                                                                                                                                       |
| Synonyms:     | MOR201-2                                                                                                                                                                                                                                                                                                            |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                             |

#### Product images:



[View online »](#)