


## Product datasheet for **GA215633**

### Mouse Olfr700 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Olfr700 CRISPRa kit - CRISPR gene activation of mouse olfactory receptor 700                                                                                                                                                                                                                                     |
| Symbol:       | Olfr700                                                                                                                                                                                                                                                                                                          |
| Components:   | GA215633G1, Olfr700 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215633G2, Olfr700 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA215633G3, Olfr700 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| Synonyms:     | MOR283-4                                                                                                                                                                                                                                                                                                         |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)