


Product datasheet for **GA215667**

Mouse Olfr235 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Olfr235 CRISPRa kit - CRISPR gene activation of mouse olfactory receptor 235
Symbol:	Olfr235
Components:	GA215667G1, Olfr235 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA215667G2, Olfr235 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA215667G3, Olfr235 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	MOR214-3; Olfr232; Olfr232-ps1; Olfr235-ps1; Olfr245; Olfr245-ps1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)