

Product datasheet for **GA215893**

Mouse **Plekhg3** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Plekhg3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology domain containing, family G (with RhoGef domain) member 3
Symbol:	Plekhg3
UniProt ID:	Q4VAC9
Synonyms:	BC030417
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)